

CharmFil Flow Blue (For base liner light curing Flowable Resin) 2x2g



CharmFil® Flow

Light Curing Flowable Composite Resin

CharmFil® Flow is a light curing flowable composite resin. It has the flexibility and elasticity which is the ideal for small cavity filling and could be applied in Cervical or Proximal. It has a high hardness and elasticity make it easy to be used in both anterior and posterior teeth area.

▲ Shade: A0, A1, A2, A3, A3.5, A4, B1, B2, C2, U0, OA2, OA3, Blue, IT, DW, A0 HV, A1 HV, A2 HV, A3 HV, A3.5 HV, B2 HV, U0 HV, Blue HV

Product Characteristics

- Excellent shade Matching with natural teeth
- Radiopaque
- Excellent workability
- Universal composite resin

Application Range

- Restoration of class I ~V cavities
- Sealing of small pit and fissure
- Undercut block-out
- Base/Liner of filling materials
- Filling of cervical

Classification

- Resin, Filling, Composite

Packaging

- CharmFil® Flow 2 g x 2 syringes, Disposable tips

Clinical procedures

1, Shade selection

2, Preparation

3, Etching

4, Cleaning and dry

5, Bonding

6, Filling

7, Light Curing

8, Polishing

CharmFil Flow Blue (For base liner light curing Flowable Resin) 2x2g

1 / 2

Rating: Not Rated Yet

[Ask a question about this product](#)